# TECNIS tórikus kalkulátor

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Johnson & Johnson VISION

Registration

Europe - Middle East - Africa English I have read and agree to the terms below.

Accept

TECNIS® intraocular lenses are currently available for sale in select markets. This application is therefore intended only for residents of those markets in which the respective TECNIS® intraocular lens(es) are available and this application may only be used for those TECNIS® intraocular lens(es) that are approved for use in your market. Using this mobile application and clicking the “Accept” button constitutes your acknowledgement of these restrictions, and your acknowledgement that you selected your proper location of residence.

TECNIS® TORIC IOL CALCULATOR LICENSE AGREEMENT, TERMS OF USE AND PRIVACY POLICY

PLEASE READ THIS USE AND LICENSE AGREEMENT (“AGREEMENT”) IN FULL

PRIVACY POLICY TERMS OF SERVICE CONTACT LRI CALCULATOR

The product-specific application is intended for the residents of a particular country or countries. As a result, the application may contain information on pharmaceuticals, medical devices and other products or uses of those products that are not approved in other countries or regions. If you are a resident of a country other than those to which the application lists, please return to surgical.jnjvision.com or contact your local Johnson & Johnson Vision affiliate to obtain the appropriate product information for your country of residence.

This application is intended for use in connection with Johnson & Johnson Vision IOLs only and should not be used in connection with any other IOLs.

©2021 Johnson & Johnson Surgical Vision, Inc.

Unless otherwise specified, all product names appearing in this program are trademarks owned by or licensed to Johnson & Johnson Vision, its subsidiaries or affiliates. No use of any Johnson & Johnson Vision trademark, trade name, or trade dress in this program may be made without the prior written authorization of Johnson & Johnson Vision, except to identify the product or service of the company.

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The screenshot displays the TECNIS Family of IOLs calculator interface. It is divided into several sections:

- Surgeon and Patient Information:** Fields for Surgeon Name, Patient Information, Eye Selection (OD/OS), and K Notation.
- Keratometry:** Fields for Surgically Induced Astigmatism (SIA), Flat K1, Steep K2, and Preop Corneal Astigmatism.
- Biometry:** Fields for Axial Length (mm), Method (Optical or Immersion), and A-constant (119.30).
- Calculator:** A list of lens options categorized by type:
  - MONOFOCAL:** TECNIS Eyhance™ Toric II DIU, TECNIS Monofocal Toric II ZCU, TECNIS Monofocal Toric ZCT.
  - PRESBYOPIA CORRECTING:** TECNIS Synergy™ Toric II DFW, TECNIS Symphony™ Toric II OptiBlue™ ZXW, TECNIS Symphony™ Toric ZXT.
  - Multifocal:** TECNIS Multifocal Toric II ZKU (+2.75 D), TECNIS Multifocal Toric II ZLU (+3.25 D).
- Visual Diagram:** A circular diagram representing the eye's axis, with markings at 90°, 135°, 45°, 180°, 225°, 270°, and 315°. The vertical axis is labeled 'N a s a l'.
- Results Table:** A table with columns for IOL Model, Orientation, Residual Astigmatism (Cylinder, Axis).

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The screenshot shows the 'TORIC CALCULATOR' interface for TECNIS Eyhance Toric II DIU lenses. The interface is organized into several functional areas:

- Surgeon and Patient Information:** Includes fields for Surgeon Name, Date (3/24/2022), Patient Information, Patient Age, Eye Selection (OD/OS), and K Notation (D/mm).
- Keratometry:** Includes fields for Surgically Induced Astigmatism (SIA), Flat K1, Steep K2, and Preop Corneal Astigmatism, along with meridian selection options.
- Biometry:** Includes Axial Length (mm), Method (Optical or Immersion), and A-constant (119.30).
- Calculation Preferences:** Includes SE IOL Power, K Index, and Refractive Cylinder Convention (Plus/Minus).
- Final Results:** A table for IOL Details and Residual Astigmatism.

Additional features include a circular diagram for meridian selection (0° to 315°), a 'Notes...' field, and navigation buttons at the bottom: Calculate Results, Clear Entries, and Print Results. The footer contains links for PRIVACY POLICY, TERMS OF SERVICE, CONTACT, and LRI CALCULATOR.

<sup>1</sup> Koch DD et al. Contribution of posterior corneal astigmatism to total corneal astigmatism. J Cataract Refract Surg. 2012 Dec;38(12):2080-7

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**TECNIS Eyhance™ Toric II DIU** TORIC CALCULATOR  
United States  
CHANGE LENS

Calculator Preferences & Support

Surgeon and Patient Information

Surgeon Name: Dr. Surgeon Date: 3/24/2022  
Patient Information: XYZ-123-ABC Patient Age: 65  
Eye Selection:  OD (Right)  OS (Left) K Notation:  D  mm

Keratometry

Surgically Induced Astigmatism (SIA): 0.35 D @ Meridian (Incision Location): 180  
Flat K1: 41.25 D Flat K1 @ Meridian: 12  
Steep K2: 43.34 D Steep K2 @ Meridian: 102  
Preop Corneal Astigmatism: 2.09 D  Include Posterior Corneal Astigmatism

Biometry

Axial Length: 24.11 mm  
Method: Optical or Immersion  
A-constant: 119.30

Calculation Preferences

SE IOL Power: 20.0 D  
K Index: 1.3375  
Refractive Cylinder Convention:  Plus  Minus

Final Results

IOL Details		Residual Astigmatism	
IOL Model	Orientation	Cylinder	Axis
<input type="radio"/> DIU300	100 °	+0.38 D	100 °
<input type="radio"/> DIU375	100 °	+0.11 D	10 °
<input type="radio"/> DIU450	100 °	+0.61 D	10 °

Calculate Results Clear Entries Print Results

PRIVACY POLICY TERMS OF SERVICE CONTACT LRI CALCULATOR

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Toric II  
DIU


Return to Calculator

Surgeon Name	Patient Information	Date	Patient Age
Dr. Surgeon	XYZ-123-ABC	3/24/2022	65

**OD (Right)**

T  
e  
m  
p  
o  
r  
a  
l

N  
a  
s  
a  
l

— (Incision Location) 180°  
— Calculated Orientation

IOL Details		Residual Refraction	
IOL Model	Orientation	Cylinder	Axis
DIU375	100 °	+0.11 D	10 °

Calculation Details		
SE IOL Power	K Index	Refractive Cylinder Convention
+20.0	1.3375	Plus

Keratometry and Biometry						
Flat K1 @ Meridian	Steep K2	SIA @ Meridian	AxL Length	Biometry Method	A-constant	PCA
41.25@12°	43.34@102°	0.35@180°	24.11	Optical or Immersion	119.30	No

IOL Details		Residual Astigmatism	
IOL Model	Orientation	Cylinder	Axis
DIU300	100 °	+0.38 D	100 °
DIU375	100 °	+0.11 D	10 °
DIU450	100 °	+0.61 D	10 °

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The screenshot displays the TECNIS Eyhance LRI Calculator interface. A white dialog box titled "Clear Entries" is overlaid on the screen, asking for confirmation to clear all entries. The dialog text reads: "Please confirm that you wish to clear all entries that you have made on this screen. (e.g., OD (Right) & OS (Left))". There are "OK" and "Cancel" buttons at the bottom of the dialog.

The background interface includes the following sections:

- Surgeon and Patient Information:** Surgeon Name (Dr. Surgeon), Date (3/24/2022), Patient Information (XYZ-123-ABC), Patient Age (65).
- Eye Selection:** OD (Right) selected, OS (Left) unselected. K Notation: D selected, mm unselected.
- Keratometry:** Surgically Induced Astigmatism (SIA) 0.35 D @ Meridian (Incision Location) 180; Flat K1 41.25 D Flat K1 @ Meridian 12; Steep K2 43.34 D Steep K2 @ Meridian 102; Preop Corneal Astigmatism 2.09 D. Includes a checkbox for "Include Posterior Corneal Astigmatism".
- Biometry:** Axial Length 24.11 mm, Method Optical or Immersion, A-constant 119.30.
- Calculation Preferences:** SE IOL Power 20.0 D, K Index 1.3375, Refractive Cylinder Convention Plus selected, Minus unselected.
- Final Results:** A table with columns for IOL Model, Orientation, Residual Astigmatism (Cylinder, Axis).

At the bottom of the interface are buttons for "Calculate Results", "Clear Entries", and "Print Results". The footer contains links for "PRIVACY POLICY", "TERMS OF SERVICE", "CONTACT", and "LRI CALCULATOR".

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**Preferences** Reset To Defaults Apply Preferences

**Display**

Location: United States | Language: English

K Notation:  D  mm

Decimal:  0.00  0,00

**Calculation**

Refractive Cylinder Convention:  Plus  Minus

K Index: -Select-

Method: Optical or Immersion

A-constant: 119.30

Include Posterior Corneal Astigmatism

**Support**

Training and Information Documents

- [TECNIS® Toric IOL Calculator User Tutorial](#)
- [TECNIS® Toric IOL Calculator FAQ](#)
- [TECNIS® Toric IOL Calculator Terms of Service](#)

(01)05050474708822(8012) - 4.5.2.0 CE

PRIVACY POLICY | TERMS OF SERVICE | CONTACT | LRI CALCULATOR

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The screenshot displays the 'Preferences & Support' section of the TECNIS TORIC CALCULATOR interface. The page header includes the TECNIS logo, 'TORIC CALCULATOR United States', and the Johnson & Johnson VISION logo. The interface is divided into two main sections: 'Preferences' and 'Support'.

**Preferences Section:**

- Display Location:** United States (dropdown menu)
- Language:** English (dropdown menu)
- K Notation:** D (selected), mm (radio button)
- Decimal:** 0.00 (selected), 0,00 (radio button)
- Calculation Refractive Cylinder Convention:** Plus (selected), Minus (radio button)
- K Index:** -Select- (dropdown menu)
- Method:** Optical or Immersion (dropdown menu)
- A-constant:** 119.30 (input field)
- Include Posterior Corneal Astigmatism

**Support Section:**

- Training and Information Documents:**
  - [TECNIS® Toric IOL Calculator User Tutorial](#)
  - [TECNIS® Toric IOL Calculator FAQ](#)
  - [TECNIS® Toric IOL Calculator Terms of Service](#)


At the bottom right of the Support section, the version number (01)05050474708822(8012) - 4.5.2.0 and the CE mark are displayed.

The footer contains four links: [PRIVACY POLICY](#), [TERMS OF SERVICE](#), [CONTACT](#), and [LRI CALCULATOR](#).

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Johnson & Johnson Surgical Vision, Inc.  
31 Technology Drive, Suite 200,  
Irvine, CA 92618 USA

EC REP

AMO Ireland  
Block B  
Liffey Valley Office Campus  
Quarryvale, Co. Dublin, Ireland

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AMO Australia Pty Ltd.  
1-5 Khartoum Road,  
North Ryde, NSW 2113  
Australia

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© 2023 Johnson & Johnson Surgical Vision, Inc.

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AMO Switzerland GmbH
A TECNIS név az Johnson & Johnson Surgical Vision, Inc. tulajdonát képező vagy általa engedéllyel használt védjegy.