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Johnson & Johnson VISION

Registration

Europe - Middle East - Africa English I have read and agree to the terms below.

Accept

TECNIS™ intraocular lenses are currently available for sale in select markets. This application is therefore intended only for residents of those markets in which the respective TECNIS™ intraocular lens(es) are available, and this application may only be used for those TECNIS™ intraocular lens(es) that are approved for use in your market. Using this mobile application and clicking the "Accept" button constitutes your acknowledgement of these restrictions, and your acknowledgement that you selected your proper location of residence.

TECNIS™ TORIC CALCULATOR LICENSE AGREEMENT, TERMS OF USE AND PRIVACY POLICY

PLEASE READ THIS USE AND LICENSE AGREEMENT ("AGREEMENT") IN FULL

TERMS OF SERVICE LRI CALCULATOR

The product-specific application is intended for the residents of a particular country or countries. As a result, the application may contain information on pharmaceuticals, medical devices and other products or uses of those products that are not approved in other countries or regions. If you are a resident of a country other than those to which the application lists, please return to surgical.jnvision.com or contact your local Johnson & Johnson Vision affiliate to obtain the appropriate product information for your country of residence.

This application is intended for use in connection with Johnson & Johnson Vision IOLs only and should not be used in connection with any other IOLs.

©2024 Johnson & Johnson Surgical Vision, Inc.

Unless otherwise specified, all product names appearing in this program are trademarks owned by or licensed to Johnson & Johnson Vision, its subsidiaries or affiliates. No use of any Johnson & Johnson Vision trademark, trade name, or trade dress in this program may be made without the prior written authorization of Johnson & Johnson Vision, except to identify the product or service of the company.

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The screenshot displays the TECNIS Toric Calculator interface. The top left corner features the TECNIS Family of IOLs logo. The main area is divided into several sections: Surgeon and Patient Information, Keratometry, and Biometry. A central panel lists available lens options, including TECNIS Eyhance™, TECNIS Toric II 1-Piece IOL, TECNIS Toric 1-Piece IOL, TECNIS PureSee™, TECNIS Synergy™, TECNIS Symphony™, and TECNIS Multifocal. A circular diagram on the right shows the orientation of the lens axis. The bottom right corner includes a table for IOL Details and Residual Astigmatism.

IOL Details		Residual Astigmatism	
IOL Model	Orientation	Cylinder	Axis

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The screenshot displays the TECNIS™ Toric Calculator interface. The top navigation bar includes the TECNIS Eyhance™ logo, the user name 'DIU', the region 'Europe - Middle East - Africa', and the Johnson & Johnson VISION logo. Below the navigation bar, there are two tabs: 'Calculator' (active) and 'Preferences & Support'. The main content area is divided into several sections:

- Surgeon and Patient Information:** Fields for Surgeon Name, Date (11/14/2024), Patient Information, and Patient Age. Eye Selection is set to OD (Right). K Notation is set to D. Units are set to mm.
- Keratometry:** Fields for Surgically Induced Astigmatism (SIA), Flat K1, Steep K2, and Preop Corneal Astigmatism. Includes checkboxes for 'Include Posterior Corneal Astigmatism'.
- Biometry:** Fields for Axial Length (mm), Method (Optical or Immersion), and A-constant (119.30).
- Calculation Preferences:** Fields for SE IOL Power, K Index, and Refractive Cylinder Convention (Plus/Minus).

On the right side, there is a diagram of a human eye with a circular grid for 'Nasal' and 'Temporal' meridians. The grid shows angles from 0° to 315° in 45° increments. Below the diagram is a 'Final Results' section with a table for IOL Details and Residual Astigmatism.

IOL Details		Residual Astigmatism	
IOL Model	Orientation	Cylinder	Axis

At the bottom of the interface, there are buttons for 'Calculate Results', 'Clear Entries', and 'Print Results'. The footer contains 'TERMS OF SERVICE' and 'LRI CALCULATOR'.

¹ Koch DD et al. Contribution of posterior corneal astigmatism to total corneal astigmatism. J Cataract Refract Surg. 2012 Dec;38(12):2080-7.

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The screenshot displays the TECNIS Toric Calculator interface. The top header includes the TECNIS Eyhance logo, the model name 'DIU', and the Johnson & Johnson VISION logo. The interface is divided into two main sections: 'Calculator' and 'Preferences & Support'.

Calculator Section:

- Surgeon and Patient Information:** Surgeon Name (Surgeon), Date (11/14/2024), Patient Information (XYZ-123-ABC), Patient Age (65). Eye Selection: OD (Right) selected, OS (Left) unselected. K Notation: D selected, mm unselected.
- Keratometry:** Surgically Induced Astigmatism (SIA) 0.30 D @ Meridian (Incision Location) 180. Flat K1 41.25 D Flat K1 @ Meridian 12. Steep K2 43.34 D Steep K2 @ Meridian 102. Preop Corneal Astigmatism 2.09 D. Include Posterior Corneal Astigmatism: unchecked.
- Biometry:** Axial Length 24.11 mm. Method: Optical or Immersion. A-constant 119.30.
- Calculation Preferences:** SE IOL Power 20.0 D. K Index 1.3375. Refractive Cylinder Convention: Plus selected, Minus unselected.

Preferences & Support Section:

- Final Results:** A circular diagram shows the incision location (blue line) and the residual astigmatism axis (orange/red line) at 101°. The diagram is labeled with 'Temporal' and 'Nasal' directions. A table below provides IOL model suggestions:

IOL Details		Residual Astigmatism	
IOL Model	Orientation	Cylinder	Axis
<input type="radio"/> DIU300	101 °	+0.33 D	101 °
<input type="radio"/> DIU375	101 °	+0.16 D	11 °
<input type="radio"/> DIU450	101 °	+0.65 D	11 °

Buttons at the bottom include 'Calculate Results', 'Clear Entries', and 'Print Results'. The footer contains 'TERMS OF SERVICE' and 'LRI CALCULATOR'.

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DIU


[Return to Calculator](#)

Surgeon Name	Patient Information	Date	Patient Age
Surgeon	XYZ-123-ABC	2/6/2024	65

OD (Right)

T
e
m
p
o
r
a
l

N
a
s
a
l

— (Incision Location) 180°

— Calculated Orientation

IOL Details		Residual Refraction	
IOL Model	Orientation	Cylinder	Axis
DIU300	101 °	+0.33 D	101 °

Calculation Details		
SE IOL Power	K Index	Refractive Cylinder Convention
+20.0	1.3375	Plus

Keratometry and Biometry						
Flat K1 @ Meridian	Steep K2	SIA @ Meridian	AxL Length	Biometry Method	A-constant	PCA
41.25@12°	43.34@102°	0.30@180°	24.11	Optical or Immersion	119.30	No

IOL Details		Residual Astigmatism	
IOL Model	Orientation	Cylinder	Axis
DIU300	101 °	+0.33 D	101 °
DIU375	101 °	+0.16 D	11 °
DIU450	101 °	+0.65 D	11 °

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The screenshot displays the TECNIS Eyhance LRI Calculator interface. A modal dialog box titled "Clear Entries" is centered on the screen, asking for confirmation to clear all data. The background interface includes the following sections:

- Surgeon and Patient:** Surgeon Name (Surgeon), Date (11/14/2024), Patient Information (XYZ-123-ABC), Patient Age (65).
- Eye Selection:** OD (Right) selected, OS (Left), K Notation (D selected, mm).
- Keratometry:** Surgically Induced Astigmatism (SIA) 0.30 D @ Meridian (Incision Location) 180; Flat K1 41.25 D Flat K1 @ Meridian 12; Steep K2 43.34 D Steep K2 @ Meridian 102; Preop Corneal Astigmatism 2.09 D. Includes checkbox for "Include Posterior Corneal Astigmatism".
- Biometry:** Axial Length 24.11 mm, Method Optical or Immersion, A-constant 119.30.
- Calculation Preferences:** SE IOL Power 20.0 D, K Index 1.3375, Refractive Cylinder Convention Plus selected, Minus.
- Final Results:** Table with columns for IOL Model, Orientation, Residual Astigmatism (Cylinder, Axis).

At the bottom, there are buttons for "Calculate Results", "Clear Entries", and "Print Results". The footer contains "TERMS OF SERVICE" and "LRI CALCULATOR".

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Include Posterior Corneal Astigmatism

Support

Training and Information Documents

- [TECNIS™ Toric Calculator User Tutorial](#)
- [TECNIS™ Toric Calculator Frequently Asked Questions \(FAQ\)](#)
- [TECNIS™ Toric Calculator Terms of Service](#)

Johnson & Johnson Surgical Vision, Inc.
31 Technology Drive, Suite 200
Irvine, CA 92618 USA
2024-11-14

AMO Ireland
Block B
Liffey Valley Office Campus
Quarryvale, Co. Dublin, Ireland

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Preferences

Display

Location: Europe - Middle East - Africa

Language: English

K Notation: D mm

Decimal: 0.00 0,00

Calculation

Refractive Cylinder Convention: Plus Minus

K Index: -Select-

Method: Optical or Immersion

A-constant: 119.30

Include Posterior Corneal Astigmatism

Support

Training and Information Documents

[TECNIS™ Toric Calculator User Tutorial](#)

[TECNIS™ Toric Calculator Frequently Asked Questions \(FAQ\)](#)

[TECNIS™ Toric Calculator Terms of Service](#)

Johnson & Johnson Surgical Vision, Inc.
31 Technology Drive, Suite 200
Irvine, CA 92618 USA
2024-11-14

AMO Ireland
Block B
Liffey Valley Office Campus
Quarryvale, Co. Dublin, Ireland

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Johnson & Johnson Surgical Vision, Inc.
31 Technology Drive, Suite 200,
Irvine, CA 92618 USA

AMO Ireland

Block B
Liffey Valley Office Campus
Quarryvale, Co. Dublin, Ireland

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AMO Australia Pty Ltd.
1-5 Khartoum Road,
North Ryde, NSW 2113
Australia

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© 2024 Johnson & Johnson Surgical Vision, Inc.

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AMO Switzerland GmbH
TECNIS è un marchio commerciale di proprietà o concesso su licenza a Johnson & Johnson Surgical Vision, Inc.