TECNIS toric calculator

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Johnson & Johnson VISION

Registration

Europe - Middle East - Africa English I have read and agree to the terms below.

Accept

TECNIS® intraocular lenses are currently available for sale in select markets. This application is therefore intended only for residents of those markets in which the respective TECNIS® intraocular lens(es) are available and this application may only be used for those TECNIS® intraocular lens(es) that are approved for use in your market. Using this mobile application and clicking the "Accept" button constitutes your acknowledgement of these restrictions, and your acknowledgement that you selected your proper location of residence.

TECNIS® TORIC IOL CALCULATOR LICENSE AGREEMENT, TERMS OF USE AND PRIVACY POLICY

PLEASE READ THIS USE AND LICENSE AGREEMENT ("AGREEMENT") IN FULL

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The product-specific application is intended for the residents of a particular country or countries. As a result, the application may contain information on pharmaceuticals, medical devices and other products or uses of those products that are not approved in other countries or regions. If you are a resident of a country other than those to which the application lists, please return to surgical.jnjvision.com or contact your local Johnson & Johnson Vision affiliate to obtain the appropriate product information for your country of residence.

This application is intended for use in connection with Johnson & Johnson Vision IOLs only and should not be used in connection with any other IOLs.

©2021 Johnson & Johnson Surgical Vision, Inc.

Unless otherwise specified, all product names appearing in this program are trademarks owned by or licensed to Johnson & Johnson Vision, its subsidiaries or affiliates. No use of any Johnson & Johnson Vision trademark, trade name, or trade dress in this program may be made without the prior written authorization of Johnson & Johnson Vision, except to identify the product or service of the company.

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TECNIS Eyhance™ Toric II DIU
United States
CHANGE LENS

Calculator Preferences & Support

Surgeon and Patient Information Notes...

Surgeon Name [] Date 3/24/2022 []
Patient Information [] Patient Age []
Eye Selection OD (Right) OS (Left) K Notation D mm

Keratometry

Surgically Induced Astigmatism (SIA) [] D @ Meridian (Incision Location) []
Flat K1 [] D Flat K1 @ Meridian []
Steep K2 [] D Steep K2 @ Meridian []
Preop Corneal Astigmatism [] Include Posterior Corneal Astigmatism

Biometry

Axial Length [] mm
Method Optical or Immersion
A-constant 119.30

Calculation Preferences

SE IOL Power -Select- D
K Index -Select-
Refractive Cylinder Convention Plus Minus

Final Results

IOL Details		Residual Astigmatism	
IOL Model	Orientation	Cylinder	Axis

Calculate Results Clear Entries Print Results

PRIVACY POLICY TERMS OF SERVICE CONTACT LRI CALCULATOR

¹ Koch DD et al. Contribution of posterior corneal astigmatism to total corneal astigmatism. J Cataract Refract Surg. 2012 Dec;38(12):2080-7

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Pagina Calculator - berekeningsresultaten

TECNIS Eyhance™ Toric II DIU TORIC CALCULATOR
United States
CHANGE LENS

Calculator Preferences & Support

Surgeon and Patient Information

Surgeon Name: Dr. Surgeon Date: 3/24/2022
Patient Information: XYZ-123-ABC Patient Age: 65
Eye Selection: OD (Right) OS (Left) K Notation: D mm

Keratometry

Surgically Induced Astigmatism (SIA): 0.35 D @ Meridian (Incision Location): 180
Flat K1: 41.25 D Flat K1 @ Meridian: 12
Steep K2: 43.34 D Steep K2 @ Meridian: 102
Preop Corneal Astigmatism: 2.09 D Include Posterior Corneal Astigmatism

Biometry

Axial Length: 24.11 mm
Method: Optical or Immersion
A-constant: 119.30

Calculation Preferences

SE IOL Power: 20.0 D
K Index: 1.3375
Refractive Cylinder Convention: Plus Minus

Calculate Results Clear Entries Print Results

Final Results

IOL Details		Residual Astigmatism	
IOL Model	Orientation	Cylinder	Axis
<input type="radio"/> DIU300	100 °	+0.38 D	100 °
<input type="radio"/> DIU375	100 °	+0.11 D	10 °
<input type="radio"/> DIU450	100 °	+0.61 D	10 °

PRIVACY POLICY TERMS OF SERVICE CONTACT LRI CALCULATOR

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Toric II
DIU


Return to Calculator

Surgeon Name	Patient Information	Date	Patient Age
Dr. Surgeon	XYZ-123-ABC	3/24/2022	65

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Pagina Calculator - Clear Entries (invoer wissen)

The screenshot displays the TECNIS Eyhance LRI Calculator interface. A white dialog box titled "Clear Entries" is centered on the screen, containing the text: "Please confirm that you wish to clear all entries that you have made on this screen. (e.g., OD (Right) & OS (Left))". Below the text are "OK" and "Cancel" buttons. The background interface is partially visible, showing various input fields and sections:

- Surgeon and Patient Information:** Surgeon Name (Dr. Surgeon), Date (3/24/2022), Patient Information (XYZ-123-ABC), Patient Age (65).
- Eye Selection:** OD (Right) selected, OS (Left) unselected. K Notation: D selected, mm unselected.
- Keratometry:** Surgically Induced Astigmatism (SIA) 0.35 D @ Meridian (Incision Location) 180. Flat K1 41.25 D Flat K1 @ Meridian 12. Steep K2 43.34 D Steep K2 @ Meridian 102. Preop Corneal Astigmatism 2.09 D. Include Posterior Corneal Astigmatism checkbox is unchecked.
- Biometry:** Axial Length 24.11 mm. Method: Optical or Immersion. A-constant 119.30.
- Calculation Preferences:** SE IOL Power 20.0 D. K Index 1.3375. Refractive Cylinder Convention: Plus selected, Minus unselected.
- Final Results:** A table with columns for IOL Model, Orientation, Residual Astigmatism (Cylinder, Axis).

At the bottom of the interface are buttons for "Calculate Results", "Clear Entries", and "Print Results". The footer contains links for "PRIVACY POLICY", "TERMS OF SERVICE", "CONTACT", and "LRI CALCULATOR".

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TECNIS® TORIC CALCULATOR
Family of IOLs United States Johnson & Johnson VISION

Calculator Preferences & Support

Preferences [Reset To Defaults](#) [Apply Preferences](#)

Display

Location: United States Language: English K Notation: D mm
Decimal: 0.00 0,00

Calculation

Refractive Cylinder Convention: Plus Minus K Index: -Select-
Method: Optical or Immersion A-constant: 118.30

Include Posterior Corneal Astigmatism

Support

Training and Information Documents
[TECNIS® Toric IOL Calculator User Tutorial](#)
[TECNIS® Toric IOL Calculator FAQ](#)
[TECNIS® Toric IOL Calculator Terms of Service](#)

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

The screenshot shows the 'Preferences & Support' page for the TECNIS TORIC CALCULATOR. The page is divided into two main sections: 'Preferences' and 'Support'. The 'Preferences' section includes settings for Display (Location, Language, K Notation, Decimal), Calculation (Refractive Cylinder Convention, K Index), and Method (Optical or Immersion, A-constant). The 'Support' section includes links to Training and Information Documents (User Tutorial, FAQ, Terms of Service) and the version number (01)05050474708822(8012) - 4.5.2.0. The page footer contains links for PRIVACY POLICY, TERMS OF SERVICE, CONTACT, and LRI CALCULATOR.

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Johnson & Johnson Surgical Vision, Inc.
31 Technology Drive, Suite 200,
Irvine, CA 92618 USA

AMO Ireland
Block B

Liffey Valley Office Campus
Quarryvale, Co. Dublin, Ireland

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AMO Australia Pty Ltd.
1-5 Khartoum Road,
North Ryde, NSW 2113
Australia

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© 2023 Johnson & Johnson Surgical Vision, Inc.

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AMO Switzerland GmbH
TECNIS is een handelsmerk dat in eigendom of licentie is bij Johnson & Johnson Surgical Vision, Inc.