# TECNIS<sup>™</sup> Toric Calculator

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## Pagina Registration (registratie)

Johnson & Johnson TECNIS™ Toric Calculator TECNIS™ Family of IOLs

Registration

Europe - Middle East - Africa English I have read and agree to the terms below. Accept

TECNIS™ intraocular lenses are currently available for sale in select markets. This application is therefore intended only for residents of those markets in which the respective TECNIS™ intraocular lens(es) are available, and this application may only be used for those TECNIS™ intraocular lens(es) that are approved for use in your market. Using this mobile application and clicking the "Accept" button constitutes your acknowledgement of these restrictions, and your acknowledgement that you selected your proper location of residence.

TECNIS™ TORIC CALCULATOR LICENSE AGREEMENT AND TERMS OF USE

PLEASE READ THIS USE AND LICENSE AGREEMENT ("AGREEMENT") IN FULL. THIS IS A LEGALLY BINDING CONTRACT BETWEEN YOU ("USER") AND JOHNSON

[TERMS OF SERVICE](#) [LRI CALCULATOR](#)

The product-specific application is intended for the residents of a particular country or countries. As a result, the application may contain information on pharmaceuticals, medical devices and other products or uses of those products that are not approved in other countries or regions. If you are a resident of a country other than those to which the application lists, please return to [surgical.jnjvision.com](http://surgical.jnjvision.com) or contact your local Johnson & Johnson Vision affiliate to obtain the appropriate product information for your country of residence.

This application is intended for use in connection with Johnson & Johnson Vision IOLs only and should not be used in connection with any other IOLs.

© Johnson & Johnson and its affiliates 2025

Unless otherwise specified, all product names appearing in this program are trademarks owned by or licensed to Johnson & Johnson Vision, its subsidiaries or affiliates. No use of any Johnson & Johnson Vision trademark, trade name, or trade dress in this program may be made without the prior written authorization of Johnson & Johnson Vision, except to identify the product or service of the company.

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## Pagina Calculator - lens selecteren

The screenshot displays the 'TECNIS™ Toric Calculator' interface. The page is titled 'Johnson & Johnson' and 'Canada'. A red bar at the top contains the word 'Calculator'. The main content area is divided into several sections:

- Surgeon and Patient Information:** Fields for Surgeon Name, Patient Information, Eye Selection (OD/OS), and K Notation.
- Keratometry:** Fields for Surgically Induced Astigmatism (SIA), Flat K1, Steep K2, and Preop Corneal Astigmatism.
- Biometry:** Fields for Axial Length, Method (Optical or Immersion), and A-constant (119.30).
- Calculator:** Fields for SE IOL Power, K Index, and Refractive Cylinder Convention (Plus/Minus).

A central dropdown menu titled 'SELECT LENS' is open, showing a list of lens options:

- TECNIS Eyhance™ DIU (with TECNIS SIMPLICITY™)
- TECNIS™ Toric II 1-Piece IOL ZCU
- TECNIS Odyssey™ DRT (with TECNIS SIMPLICITY™)
- TECNIS PureSee™ DET (with TECNIS SIMPLICITY™)
- TECNIS Synergy™ DFW (with TECNIS SIMPLICITY™, OptiBlue™)
- TECNIS Symphony™ ZXT

On the right side, there is a 'Preferences & Support' section with a diagram of an eye and a circular scale for axis selection (0°, 45°, 90°, 135°, 180°, 225°, 270°, 315°). Below this is a 'Results' section with a table for 'IOL Details' and 'Residual Astigmatism'.

At the bottom, there are buttons for 'Calculate Results', 'Clear Entries', and 'Print Results'. A red bar at the very bottom contains 'TERMS OF SERVICE' and 'LRI CALCULATOR'.

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The screenshot displays the TECNIS™ Toric Calculator interface. The top header includes the Johnson & Johnson logo, the product name 'TECNIS™ Toric Calculator', the region 'Europe - Middle East - Africa', and the 'CHANGE LENS' button. The right side features the 'TECNIS Eyhance™ with TECNIS SIMPLICITY™' logo and 'DIU'.

The main interface is divided into two tabs: 'Calculator' (active) and 'Preferences & Support'. The 'Calculator' tab contains several sections:

- Surgeon and Patient Information:** Includes fields for Surgeon Name, Date (7/29/2025), Patient Information, Patient Age, Eye Selection (OD (Right) selected), K Notation (D selected), and units (mm selected).
- Keratometry:** Includes fields for Surgically Induced Astigmatism (SIA), Flat K1, Steep K2, and Preop Corneal Astigmatism. It also has checkboxes for 'Include Posterior Corneal Astigmatism' and 'Include Meridian (Incision Location)'. There are also fields for 'Flat K1 @ Meridian' and 'Steep K2 @ Meridian'.
- Biometry:** Includes fields for Axial Length (mm), Method (Optical or Immersion), and A-constant (119.30).
- Calculation Preferences:** Includes dropdowns for SE IOL Power and K Index, and radio buttons for Refractive Cylinder Convention (Plus selected, Minus).

At the bottom of the calculator section are buttons for 'Calculate Results', 'Clear Entries', and 'Print Results'. To the right of the input fields is a diagram of a human eye with a circular meridian scale. The scale is labeled 'T e m p o r a l' on the left and 'N a s a l' on the right. The scale has markings at 90°, 45°, 0°, 315°, 270°, 225°, and 180°. A blue box highlights the 'OD' (Right) eye selection.

The 'Final Results' section contains a table with the following structure:

IOL Details		Residual Astigmatism	
IOL Model	Orientation	Cylinder	Axis

At the bottom of the page are links for 'TERMS OF SERVICE' and 'LRI CALCULATOR'.

<sup>1</sup> Koch DD et al. Contribution of posterior corneal astigmatism to total corneal astigmatism. J Cataract Refract Surg. 2012 Dec;38(12):2080-7

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## Pagina Calculator - berekeningsresultaten

The screenshot displays the TECNIS Toric Calculator interface. The top navigation bar includes the Johnson & Johnson logo, the product name 'TECNIS™ Toric Calculator', regional information 'Europe - Middle East - Africa', and the IOL model 'TECNIS Eyhance™ with TECNIS SIMPLICITY™ DIU'. A 'CHANGE LENS' button is visible.

The interface is divided into two main sections: 'Calculator' and 'Preferences & Support'. The 'Calculator' section contains several input panels:

- Surgeon and Patient Information:** Surgeon Name (Surgeon), Date (7/29/2025), Patient Information (XYZ-123-ABC), Patient Age (65), Eye Selection (OD (Right)), K Notation (D), and units (mm).
- Keratometry:** Surgically Induced Astigmatism (SIA) (0.30 D), Flat K1 (41.25 D), Steep K2 (43.34 D), Preop Corneal Astigmatism (2.09 D), and Meridian Location (180°).
- Biometry:** Axial Length (24.11 mm), Method (Optical or Immersion), and A-constant (119.30).
- Calculation Preferences:** SE IOL Power (20.0 D), K Index (1.3375), and Refractive Cylinder Convention (Plus).

The 'Preferences & Support' section features a diagram of an eye with a blue box highlighting the incision location and an orange/red line indicating the astigmatism axis. The diagram is labeled 'OD' and includes 'Temporal' and 'Nasal' markers. The angle of the astigmatism axis is shown as 101°.

The 'Final Results' section displays a table with the following data:

IOL Details		Residual Astigmatism	
IOL Model	Orientation	Cylinder	Axis
DIU300	101°	+0.33 D	101°
DIU375	101°	+0.16 D	11°
DIU450	101°	+0.65 D	11°

Buttons for 'Calculate Results', 'Clear Entries', and 'Print Results' are located at the bottom of the calculator section. The footer contains 'TERMS OF SERVICE' and 'LRI CALCULATOR' links.

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## Pagina Results (resultaten) afdrukken

[Return to Calculator](#)


**TECNIS**  
**Eyhance™**  
with TECNIS SIMPLICITY™

DIU

Surgeon Name	Patient Information	Date	Patient Age
Surgeon	XYZ-123-ABC	7/30/2025	65

**OD (Right)**

T  
e  
m  
p  
o  
r  
a  
l

N  
a  
s  
a  
l

— (Incision Location) 180°

— Calculated Orientation

IOL Details		Residual Refraction	
<b>IOL Model</b>	<b>Orientation</b>	<b>Cylinder</b>	<b>Axis</b>
DIU300	101 °	+0.33 D	101 °

Calculation Details		
<b>SE IOL Power</b>	<b>K Index</b>	<b>Refractive Cylinder Convention</b>
+20.0	1.3375	Plus

Keratometry and Biometry						
<b>Flat K1 @ Meridian</b>	<b>Steep K2</b>	<b>SIA @ Meridian</b>	<b>AxL Length</b>	<b>Biometry Method</b>	<b>A-constant</b>	<b>PCA</b>
41.25@12°	43.34@102°	0.30@180°	24.11	Optical or Immersion	119.30	No

IOL Details		Residual Astigmatism	
<b>IOL Model</b>	<b>Orientation</b>	<b>Cylinder</b>	<b>Axis</b>
DIU300	101 °	+0.33 D	101 °
DIU375	101 °	+0.16 D	11 °
DIU450	101 °	+0.65 D	11 °

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## Pagina Calculator - Clear Entries (invoer wissen)

The screenshot displays the Johnson & Johnson LRI Calculator interface. A modal dialog box titled "Clear Entries" is centered on the screen, asking for confirmation to clear all entries. The background interface includes the following sections:

- Surgeon and Patient Information:** Surgeon Name, Date (7/29/2025), Patient Information (XYZ-123-ABC), Patient Age (65).
- Eye Selection:** OD (Right) selected, OS (Left), K Notation (D selected, mm).
- Keratometry:** Surgically Induced Astigmatism (SIA) 0.30 D @ Meridian (Incision Location) 180; Flat K1 41.25 D Flat K1 @ Meridian 12; Steep K2 43.34 D Steep K2 @ Meridian 102; Preop Corneal Astigmatism 2.09 D. Includes checkbox for "Include Posterior Corneal Astigmatism".
- Biometry:** Axial Length 24.11 mm, Method Optical or Immersion, A-constant 119.30.
- Calculation Preferences:** SE IOL Power 20.0 D, K Index 1.3375, Refractive Cylinder Convention Plus selected, Minus.
- Final Results:** Table with columns for IOL Model, Orientation, Residual Astigmatism (Cylinder, Axis).

At the bottom of the interface are buttons for "Calculate Results", "Clear Entries", and "Print Results". The footer contains "TERMS OF SERVICE" and "LRI CALCULATOR".

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The screenshot shows the 'Preferences & Support' page for the TECNIS™ Toric Calculator. The page is divided into two main sections: 'Preferences' and 'Support'. The 'Preferences' section is highlighted with a blue border and contains the following settings:

- Display:**
  - Location: Europe - Middle East - Africa
  - Language: English
  - K Notation: D (selected), mm
  - Decimal: 0.00 (selected), 0,00
- Calculation:**
  - Refractive Cylinder Convention: Plus (selected), Minus
  - K Index: -Select-
  - Method: Optical or Immersion
  - A-constant: 119.30
  - Include Posterior Corneal Astigmatism

The 'Support' section includes links for 'Training and Information Documents', 'TECNIS™ Toric Calculator User Tutorial', 'TECNIS™ Toric Calculator Frequently Asked Questions (FAQ)', and 'TECNIS™ Toric Calculator Terms of Service'. It also provides contact information for Johnson & Johnson Surgical Vision, Inc. and regulatory information.

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The screenshot shows the 'Preferences & Support' page for the TECNIS™ Toric Calculator. The page is divided into two main sections: 'Preferences' and 'Support'. The 'Preferences' section is currently active and shows settings for 'Display' (Location: Europe - Middle East - Africa, Language: English, K Notation: mm, Decimal: 0.00) and 'Calculation' (Refractive Cylinder Convention: Plus, K Index: -Select-, A-constant: 119.30). The 'Support' section is highlighted in blue and contains links for 'Training and Information Documents', 'TECNIS™ Toric Calculator User Tutorial', 'TECNIS™ Toric Calculator Frequently Asked Questions (FAQ)', and 'TECNIS™ Toric Calculator Terms of Service'. It also includes contact information for Johnson & Johnson Surgical Vision, Inc. and AMO Ireland, along with regulatory information (MD, CE, LUDI) and a phone number.

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Johnson & Johnson Surgical Vision, Inc.  
31 Technology Drive, Suite 200,  
Irvine, CA 92618 USA

AMO Ireland  
Block B

Liffey Valley Office Campus  
Quarryvale, Co. Dublin, Ireland

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AMO Australia Pty Ltd.  
1-5 Khartoum Road,  
North Ryde, NSW 2113  
Australia

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AMO Switzerland GmbH
TECNIS is een handelsmerk dat in eigendom of licentie is bij Johnson & Johnson Surgical Vision, Inc.