Calculador tórico TECNIS

Manual de funcionamento

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The screenshot displays the TECNIS Family of IOLs calculator interface. The main content area is divided into three sections: MONOFOCAL, PRESBYOPIA CORRECTING, and a list of selected lenses. The MONOFOCAL section lists TECNIS Eyhance™ DIU (Toric II), TECNIS Monofocal ZCU (Toric II), and TECNIS Monofocal ZCT (Toric). The PRESBYOPIA CORRECTING section lists TECNIS Synergy™ DFW (Toric II), TECNIS Symphony™ ZXW (Toric II, OptiBlue™), TECNIS Symphony™ ZXT (Toric), TECNIS Multifocal ZKU (Toric II, +2.75 D), and TECNIS Multifocal ZLU (Toric II, +3.25 D). The interface includes input fields for Surgeon and Patient Information, Keratometry (Surgically Induced Astigmatism, Flat K1, Steep K2, Preop Corneal Astigmatism), and Biometry (Axial Length, Method, A-constant). A circular diagram on the right shows meridian angles (90°, 45°, 0°, 315°, 270°, 225°, 180°, 135°) and the text 'N a s a i'. The bottom of the interface features a red bar with links for PRIVACY POLICY, TERMS, CONTACT, and LRI CALCULATOR.

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The screenshot displays the 'TORIC CALCULATOR' interface for TECNIS Eyhance. The top navigation bar includes the TECNIS Eyhance logo, 'Toric II DIU' label, 'United States' location, and a 'CHANGE LENS' dropdown. The main interface is split into two tabs: 'Calculator' (active) and 'Preferences & Support'. The 'Calculator' tab contains several input sections: 'Surgeon and Patient Information' (with a 'Notes...' button), 'Keratometry', 'Biometry', and 'Calculation Preferences'. The 'Final Results' section is partially visible on the right, showing a table with columns for 'IOL Model', 'Orientation', 'Cylinder', and 'Axis'. To the right of the input fields is a circular diagram representing meridians (0° to 315°) and a stylized eye icon with 'OD' (Right) and 'OS' (Left) labels. The bottom of the page features a red bar with links for 'PRIVACY POLICY', 'TERMS OF SERVICE', 'CONTACT', and 'LRI CALCULATOR'.

¹ Koch DD et al. Contribution of posterior corneal astigmatism to total corneal astigmatism. J Cataract Refract Surg. 2012 Dec;38(12):2080-7

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The screenshot displays the TECNIS Eyhance TORIC CALCULATOR interface. The top navigation bar includes the TECNIS Eyhance logo, 'Toric II DIU', 'United States', and 'CHANGE LENS'. The main interface is divided into several sections:

- Surgeon and Patient Information:** Surgeon Name (Dr. Surgeon), Date (3/24/2022), Patient Information (XYZ-123-ABC), Patient Age (65), Eye Selection (OD (Right)), K Notation (D), and mm.
- Keratometry:** Surgically Induced Astigmatism (SIA) (0.35 D), Flat K1 (41.25 D), Steep K2 (43.34 D), Preop Corneal Astigmatism (2.09 D), Meridian Location (180), Flat K1 @ Meridian (12), Steep K2 @ Meridian (102), and Include Posterior Corneal Astigmatism (checkbox).
- Biometry:** Axial Length (24.11 mm), Method (Optical or Immersion), and A-constant (119.30).
- Calculation Preferences:** SE IOL Power (20.0 D), K Index (1.3375), and Refractive Cylinder Convention (Plus).
- Final Results:** A table showing IOL Model, Orientation, Residual Astigmatism (Cylinder), and Axis.

The 'Final Results' table is highlighted with a blue border:

IOL Details		Residual Astigmatism	
IOL Model	Orientation	Cylinder	Axis
<input type="radio"/> DIU300	100 °	+0.38 D	100 °
<input type="radio"/> DIU375	100 °	+0.11 D	10 °
<input type="radio"/> DIU450	100 °	+0.61 D	10 °

The interface also features a diagram of an eye with a circular scale for orientation, showing a 100° angle and various meridian markings (90°, 45°, 0°, 315°, 270°, 225°, 180°, 135°). The diagram is labeled 'OD' and 'Nasal'.

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Toric II
DIU


Return to Calculator

Surgeon Name	Patient Information	Date	Patient Age
Dr. Surgeon	XYZ-123-ABC	3/24/2022	65

OD (Right)

T
e
m
p
o
r
a
l

N
a
s
a
l

— (Incision Location) 180°
— Calculated Orientation

IOL Details		Residual Refraction	
IOL Model	Orientation	Cylinder	Axis
DIU375	100 °	+0.11 D	10 °

Calculation Details		
SE IOL Power	K Index	Refractive Cylinder Convention
+20.0	1.3375	Plus

Keratometry and Biometry						
Flat K1 @ Meridian	Steep K2	SIA @ Meridian	AxL Length	Biometry Method	A-constant	PCA
41.25@12°	43.34@102°	0.35@180°	24.11	Optical or Immersion	119.30	No

IOL Details		Residual Astigmatism	
IOL Model	Orientation	Cylinder	Axis
DIU300	100 °	+0.38 D	100 °
DIU375	100 °	+0.11 D	10 °
DIU450	100 °	+0.61 D	10 °

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The screenshot displays the TECNIS Eyhance LRI Calculator interface. A white dialog box titled "Clear Entries" is centered on the screen, containing the text: "Please confirm that you wish to clear all entries that you have made on this screen. (e.g., OD (Right) & OS (Left))". The dialog has "OK" and "Cancel" buttons. The background interface is dimmed and includes the following sections:

- Surgeon and Patient Information:** Surgeon Name (Dr. Surgeon), Date (3/24/2022), Patient Information (XYZ-123-ABC), Patient Age (65), Eye Selection (OD (Right) selected), K Notation (D selected).
- Keratometry:** Surgically Induced Astigmatism (SIA) (0.35 D), Flat K1 (41.25 D), Steep K2 (43.34 D), Preop Corneal Astigmatism (2.09 D). Includes Meridian (Incision Location) and Meridian values.
- Biometry:** Axial Length (24.11 mm), Method (Optical or Immersion), A-constant (119.30).
- Calculation Preferences:** SE IOL Power (20.0 D), K Index (1.3375), Refractive Cylinder Convention (Plus selected).
- Final Results:** A table for IOL Details and Residual Astigmatism.

At the bottom of the interface are buttons for "Calculate Results", "Clear Entries", and "Print Results". The footer contains links for "PRIVACY POLICY", "TERMS OF SERVICE", "CONTACT", and "LRI CALCULATOR".

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The screenshot shows the 'Preferences & Support' page of the TECNIS TORIC CALCULATOR. The page is divided into two main sections: 'Preferences' and 'Support'. The 'Preferences' section is further divided into 'Display', 'Calculation', and 'Method'.

Display:

- Location: United States (dropdown)
- Language: English (dropdown)
- K Notation: D, mm
- Decimal: 0.00, 0,00

Calculation:

- Refractive Cylinder Convention: Plus, Minus
- K Index: -Select- (dropdown)

Method:

- Method: Optical or Immersion (dropdown)
- A-constant: 119.30 (input field)

Include Posterior Corneal Astigmatism

Support:

- Training and Information Documents
 - [TECNIS® Toric IOL Calculator User Tutorial](#)
 - [TECNIS® Toric IOL Calculator FAQ](#)
 - [TECNIS® Toric IOL Calculator Terms of Service](#)

Footer: (01)05050474708822(8012) - 4.5.2.0 ©

Bottom navigation: PRIVACY POLICY, TERMS OF SERVICE, CONTACT, LRI CALCULATOR

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The screenshot displays the 'Preferências e suporte' page of the TECNIS TORIC CALCULATOR. The page is divided into two main sections: 'Preferences' and 'Support'. The 'Preferences' section includes settings for Display (Location: United States, Language: English), K Notation (D, mm), Decimal (0.00, 0,00), Calculation (Refractive Cylinder Convention: Plus, Minus), K Index (-Select-), Method (Optical or Immersion), and A-constant (119.30). There is also a checkbox for 'Include Posterior Corneal Astigmatism'. The 'Support' section lists 'Training and Information Documents' with links to the User Tutorial, FAQ, and Terms of Service. The footer contains the version number (01)05050474708822(8012) - 4.5.2.0 and a CE mark. The bottom navigation bar includes links for PRIVACY POLICY, TERMS OF SERVICE, CONTACT, and LRI CALCULATOR.

TECNIS®
Family of IOLs

TORIC CALCULATOR
United States

Johnson & Johnson VISION

Calculator Preferences & Support

Preferences Reset To Defaults Apply Preferences

Display
Location: United States Language: English
K Notation: D mm
Decimal: 0.00 0,00

Calculation
Refractive Cylinder Convention: Plus Minus
K Index: -Select-
Method: Optical or Immersion
A-constant: 119.30
 Include Posterior Corneal Astigmatism

Support

Training and Information Documents
[TECNIS® Toric IOL Calculator User Tutorial](#)
[TECNIS® Toric IOL Calculator FAQ](#)
[TECNIS® Toric IOL Calculator Terms of Service](#)

(01)05050474708822(8012) - 4.5.2.0 CE

PRIVACY POLICY TERMS OF SERVICE CONTACT LRI CALCULATOR

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Johnson & Johnson Surgical Vision, Inc.
31 Technology Drive, Suite 200,
Irvine, CA 92618 USA

AMO Ireland
Block B
Liffey Valley Office Campus
Quarryvale, Co. Dublin, Ireland

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AMO Australia Pty Ltd.
1-5 Khartoum Road,
North Ryde, NSW 2113
Australia

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© 2023 Johnson & Johnson Surgical Vision, Inc.

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AMO Switzerland GmbH
TECNIS é uma marca comercial pertencente ou licenciada à Johnson & Johnson Surgical Vision, Inc.