TECNIS Torik Hesaplayıcı

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TECNIS® TORIC CALCULATOR *Johnson & Johnson VISION*
Family of IOLs

Registration

Europe - Middle East - Africa English I have read and agree to the terms below.
Accept

TECNIS® intraocular lenses are currently available for sale in select markets. This application is therefore intended only for residents of those markets in which the respective TECNIS® intraocular lens(es) are available and this application may only be used for those TECNIS® intraocular lens(es) that are approved for use in your market. Using this mobile application and clicking the "Accept" button constitutes your acknowledgement of these restrictions, and your acknowledgement that you selected your proper location of residence.

TECNIS® TORIC IOL CALCULATOR LICENSE AGREEMENT, TERMS OF USE AND PRIVACY POLICY

PLEASE READ THIS USE AND LICENSE AGREEMENT ("AGREEMENT") IN FULL

PRIVACY POLICY TERMS OF SERVICE CONTACT LRI CALCULATOR

The product-specific application is intended for the residents of a particular country or countries. As a result, the application may contain information on pharmaceuticals, medical devices and other products or uses of those products that are not approved in other countries or regions. If you are a resident of a country other than those to which the application lists, please return to surgical.jnjvision.com or contact your local Johnson & Johnson Vision affiliate to obtain the appropriate product information for your country of residence.

This application is intended for use in connection with Johnson & Johnson Vision IOLs only and should not be used in connection with any other IOLs.

©2021 Johnson & Johnson Surgical Vision, Inc.

Unless otherwise specified, all product names appearing in this program are trademarks owned by or licensed to Johnson & Johnson Vision, its subsidiaries or affiliates. No use of any Johnson & Johnson Vision trademark, trade name, or trade dress in this program may be made without the prior written authorization of Johnson & Johnson Vision, except to identify the product or service of the company.

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The screenshot displays the TECNIS IOL Calculator interface. The top left features the TECNIS Family of IOLs logo. The main area is divided into sections for patient information, keratometry, and biometry. The lens selection panel on the right lists various options under MONOFOCAL and PRESBYOPIA CORRECTING categories. The bottom right shows a circular diagram for astigmatism correction with angles from 0° to 315°.

TECNIS® Family of IOLs

Calculator

Surgeon and Patient Information

Surgeon Name: [Input Field] Date: [Input Field]

Patient Information: [Input Field] Pat: [Input Field]

Eye Selection: OD (Right) OS (Left) K Notation: [Input Field]

Keratometry

Surgically Induced Astigmatism (SIA): [Input Field] D @ Meridian Location)

Flat K1: [Input Field] D Flat K1 @ M

Steep K2: [Input Field] D Steep K2 @

Preop Corneal Astigmatism: [Input Field] D Include

Biometry

Axial Length: [Input Field] mm

Method: Optical or Immersion

A-constant: 119.30

Calculate Results **Clear Entries**

MONOFOCAL

TECNIS Eyhance™ Toric II DIU
with TECNIS SIMPLICITY™ Delivery System

TECNIS Monofocal Toric II ZCU

TECNIS Monofocal Toric ZCT

PRESBYOPIA CORRECTING

TECNIS Synergy™ Toric II DFW
with TECNIS SIMPLICITY™ Delivery System

TECNIS Symphony™ Toric II OptiBlue® ZXW
Extended Range of Vision IOL

TECNIS Symphony™ Toric ZXT
Extended Range of Vision IOL

TECNIS Multifocal Toric II ZKU
+2.75 D

TECNIS Multifocal Toric II ZLU
+3.25 D

Johnson & Johnson VISION

Preferences & Support

Diagram: A circular diagram representing astigmatism correction with angles: 0°, 45°, 90°, 135°, 180°, 225°, 270°, 315°.

Results

IOL Details		Residual Astigmatism	
IOL Model	Orientation	Cylinder	Axis

PRIVACY POLICY **TERMS** **ACT** **LRI CALCULATOR**

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TECNIS Eyhance™ Toric II DIU
with TECNIS SIMPLICITY™ Delivery System

TORIC CALCULATOR
United States
CHANGE LENS

Johnson & Johnson VISION

Calculator

Preferences & Support

Surgeon and Patient Information

Surgeon Name: Dr. Surgeon Date: 3/24/2022
Patient Information: XYZ-123-ABC Patient Age: 65
Eye Selection: OD (Right) OS (Left) K Notation: D mm

Keratometry

Surgically Induced Astigmatism (SIA): 0.35 D @ Meridian (Incision Location): 180
Flat K1: 41.25 D Flat K1 @ Meridian: 12
Steep K2: 43.34 D Steep K2 @ Meridian: 102
Preop Corneal Astigmatism: 2.09 D Include Posterior Corneal Astigmatism

Biometry

Axial Length: 24.11 mm
Method: Optical or Immersion
A-constant: 119.30

Calculation Preferences

SE IOL Power: 20.0 D
K Index: 1.3375
Refractive Cylinder Convention: Plus Minus

Calculate Results Clear Entries Print Results

Final Results

IOL Details		Residual Astigmatism	
IOL Model	Orientation	Cylinder	Axis
<input type="radio"/> DIU300	100 °	+0.38 D	100 °
<input type="radio"/> DIU375	100 °	+0.11 D	10 °
<input type="radio"/> DIU450	100 °	+0.61 D	10 °

PRIVACY POLICY TERMS OF SERVICE CONTACT LRI CALCULATOR

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TECNIS Eyhance™ Toric II DIU with TECNIS SIMPLICITY™ Delivery System

Return to Calculator

Johnson & Johnson VISION

Surgeon Name	Patient Information	Date	Patient Age
Dr. Surgeon	XYZ-123-ABC	3/24/2022	65

OD (Right)

IOL Details		Residual Refraction	
IOL Model	Orientation	Cylinder	Axis
DIU375	100 °	+0.11 D	10 °

Calculation Details		
SE IOL Power	K Index	Refractive Cylinder Convention
+20.0	1.3375	Plus

Keratometry and Biometry						
Flat K1 @ Meridian	Steep K2	SIA @ Meridian	AxL Length	Biometry Method	A-constant	PCA
41.25@12°	43.34@102°	0.35@180°	24.11	Optical or Immersion	119.30	No

IOL Details		Residual Astigmatism	
IOL Model	Orientation	Cylinder	Axis
DIU300	100 °	+0.38 D	100 °
DIU375	100 °	+0.11 D	10 °
DIU450	100 °	+0.61 D	10 °

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The screenshot displays the TECNIS Eyhance IOL Calculator interface. A modal dialog box titled "Clear Entries" is open, asking for confirmation to clear all entries. The dialog text reads: "Please confirm that you wish to clear all entries that you have made on this screen. (e.g., OD (Right) & OS (Left))". There are "OK" and "Cancel" buttons. The background interface includes the following sections:

- Surgeon and Patient Information:** Surgeon Name (Dr. Surgeon), Date (3/24/2022), Patient Information (XYZ-123-ABC), Patient Age (65).
- Eye Selection:** OD (Right) selected, OS (Left) unselected, K Notation (D) selected, mm unselected.
- Keratometry:** Surgically Induced Astigmatism (SIA) (0.35 D @ Meridian (Incision Location) 180), Flat K1 (41.25 D Flat K1 @ Meridian 12), Steep K2 (43.34 D Steep K2 @ Meridian 102), Preop Corneal Astigmatism (2.09 D). Includes a checkbox for "Include Posterior Corneal Astigmatism".
- Biometry:** Axial Length (24.11 mm), Method (Optical or Immersion), A-constant (119.30).
- Calculation Preferences:** SE IOL Power (20.0 D), K Index (1.3375), Refractive Cylinder Convention (Plus selected, Minus unselected).
- Final Results:** IOL Details (IOL Model, Orientation) and Residual Astigmatism (Cylinder, Axis).

At the bottom of the interface, there are buttons for "Calculate Results", "Clear Entries", and "Print Results". The footer contains links for "PRIVACY POLICY", "TERMS OF SERVICE", "CONTACT", and "LRI CALCULATOR".

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Preferences Reset To Defaults Apply Preferences

Display

Location: United States Language: English K Notation: D mm
Decimal: 0.00 0,00

Calculation

Refractive Cylinder Convention: Plus Minus K Index: -Select-
Method: Optical or Immersion A-constant: 119.30

Include Posterior Corneal Astigmatism

Support

Training and Information Documents
[TECNIS® Toric IOL Calculator User Tutorial](#)
[TECNIS® Toric IOL Calculator FAQ](#)
[TECNIS® Toric IOL Calculator Terms of Service](#)

(01)05050474708822(8012) - 4.5.2.0 CE

PRIVACY POLICY TERMS OF SERVICE CONTACT LRI CALCULATOR

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Preferences Reset To Defaults Apply Preferences

Display

Location: United States Language: English K Notation: D mm
Decimal: 0.00 0,00

Calculation

Refractive Cylinder Convention: Plus Minus K Index: -Select-
Method: Optical or Immersion A-constant: 119.30
 Include Posterior Corneal Astigmatism

Support

Training and Information Documents

[TECNIS® Toric IOL Calculator User Tutorial](#)
[TECNIS® Toric IOL Calculator FAQ](#)
[TECNIS® Toric IOL Calculator Terms of Service](#)


(01)05050474708822(8012) - 4.5.2.0 CE

PRIVACY POLICY TERMS OF SERVICE CONTACT LRI CALCULATOR

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EC REP

CH REP

Johnson & Johnson Surgical Vision, Inc.
31 Technology Drive, Suite 200,
Irvine, CA 92618 USA

AMO Ireland
Block B
Liffey Valley Office Campus
Quarryvale, Co. Dublin, Ireland

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AMO Australia Pty Ltd.
1-5 Khartoum Road,
North Ryde, NSW 2113
Australia

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© 2023 Johnson & Johnson Surgical Vision, Inc.

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AMO Switzerland GmbH
TECNIS; Johnson & Johnson Surgical Vision, Inc. şirketine ait veya lisanslı bir ticari markadır.